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The purpose of this test is to evaluate if the tested material is or is not a skin irritant, according to the method described in OECD 439. This method uses human artificial skin model (EPISKIN™, EpiDerm™, SkinEthic™ RHE, epiCS®) to assess the skin irritation of chemical substances, mixture and to properly label them to this respect if applicable.

The test is based on the evaluation of cell survival after the exposure to the substance through MTT assay and by comparison with epidermis treated with phosphate buffer only (negative control). The MTT method is a colorimetric assay that allows to determinate the percentage of cells alive within an in vitro cultured tissue. This assay is based on the ability of the mitochondrial succinate dehydrogenase enzyme to metabolise the nitro-blue tetrazolium salt, giving a coloured compound that can be measured by spectrophotometer reading.

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The test has been carried out according to the method described in OECD 439 and in the MatTek validated protocol (1-2).

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
A reconstructed artificial human skin model comprising normal human epidermal keratinocytes, growing as an integrated three-dimensional cell culture model, perfectly mimicking the human skin *in vitro*. The model exhibits normal barrier functions (presence of a well-differentiated *stratum corneum*). It was supplied by MatTek (Bratislava, Batch 34123).

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The sample is a solid that has been tested by cutting out small disks the size of the epidermis; these were then applied to the epidermis previously wet with 25 µl of water to improve contact.

The positive control SLS is dissolved at 5% in water. PBS alone has been used as negative control.

30µl of controls have been applied on each epidermis units.

Sample and controls have been applied on each epidermis unit in three replicates. The exposure has been carried out for 35' at 37°C and 25' at room temperature. At the end of the exposure period the product is removed with multiple washings with PBS and the tissues were further incubated at 37°C, 5% CO₂ for 42 h. After the first 24 h, the medium has been discharged and renewed. At the end of the exposure, the viability assay is performed to evaluate the cell survival in the epidermis units.

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Epidermis units are treated with 1mg/ml MTT solution (3-[4,5-dimethylthiazol-2-yl]-2,5-diphenyl tetrazolium bromide) for 3 h at 37°C. The solution is then removed and replaced with isopropanol, with further 2 h incubation at room temperature. 2 aliquots of every sample are transferred to a 96 well plate for the reading. The absorbance is read at the wavelength of 570 nm with a colorimeter (Tecan Infinite 200 PRO model) equipped with a microplate reader. The absorbance values are corrected by subtracting the reading of the blanks (diluent only).

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The results are expressed in terms of viability:

% of cell viability = [mean OD(570 nm) test product / mean OD(570 nm) negative control] x 100

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for negative control (CN): the mean OD value of the 3 tissues has to be ≥ 0.8 e ≤ 2.8
 the standard deviation has to be $\leq 18\%$

for positive control (CP): the viability mean (expressed as % of the NC) has to be $\leq 20\%$
 the standard deviation has to be $< 18\%$

for the sample : the standard deviation has to be $< 18\%$

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www.abich.it

4 **PARTE QUARTA/PART FOUR** **RIFERIMENTI/REFERENCES**

1. OECD (2020), *Test No. 439: In Vitro Skin Irritation: Reconstructed Human Epidermis Test Method*, OECD Guidelines for the Testing of Chemicals, Section 4, OECD Publishing, Paris
2. In Vitro Epiderm™ Skin Irritation Test (EPI-200-SIT) For use with MatTek Corporation's Reconstructed Human Epidermal Model Epiderm. MatTek Corporation Protocol 10/02/2019.

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MatTek Corporation

EpiDerm QC (EPI-200)

LOT 34123
TESTED Post Refrigerated Storage
COMMENTS NO

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www.abich.it

Certificate of Analysis

Product: EpiDerm™ Reconstructed Human Epidermis

Lot Number: **34123**

Part#: EPI-200, EPI-212, EPI-218

Description: Reconstructed human epidermis tissue containing normal human keratinocytes. This product is for research use only. Not for use in animals, humans or diagnostic purposes.

I. Cell source

All cells used to produce EpiDerm™ are purchased or derived from tissue obtained by MatTek Corporation from accredited institutions. In all cases, consent was obtained by these institutions from the donor or the donor's legal next of kin, for use of the tissues or derivatives of the tissue for research purposes.

Keratinocyte Strain: **00267**

II. Analysis for potential biological contaminants

The cells used to produce EpiDerm™ tissue are screened for potential biological contaminants. Tests for each potential biological contaminant listed below were performed according to the test method given. Results of "Not detected" indicate that testing for the potential biological contaminant was not observed as determined by the stated test method.

HIV-1 virus – Oligonucleotide-directed amplification	Not detected
Hepatitis B virus – Oligonucleotide- directed amplification	Not detected
Hepatitis C virus – Oligonucleotide- directed amplification	Not detected
Bacteria, yeast, and other fungi – long term antibiotic, antimycotic free culture	Not detected

III. Analysis for tissue functionality and quality

Test	Specification	Acceptance criteria	Result and QA Statement	
Tissue viability	MTT QC assay, 4 hours, n=3	OD (540-570 nm) [1.0-3.0]	1.813 ± 0.11	Pass
Barrier function	ET-50 assay, 100 µL 1% Triton X-100, 4 time-points, n=3, MTT assay	ET-50 [4.77-8.72 hrs]	6.59 hrs	Pass
Sterility	Long term antibiotic and antimycotic free culture	No contamination	Sterile	Pass

Tissue viability and the barrier function test are within the acceptable ranges and indicate appropriate formation of the epidermal barrier, the presence of a functional stratum corneum, a viable basal cell layer, and intermediate spinous and granular layers. Results obtained with this lot conform to the requirements of the **OECD TG 431 and 439**.

Initials: **SL**
 Date: **17.02.2021**

Nelson Rivas
 Quality Assurance Associate

February 17, 2021
 Date

CAUTION: Whereas all information herein is believed to be correct, no absolute guarantee that human derived material is non-infectious can be made or is implied by this certificate of analysis. All tissues should be treated as potential pathogens. The use of protective clothing and eyewear and appropriate disposal procedures are strongly recommended.

MatTek In Vitro Life Science Laboratories
 Mlynské Nivy 73, Bratislava - Slovakia
 +421-2-3260-7401

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